

MANUAL

CUSTOMIZE
ACTIVITIES



Activity types

Games



Generators




Worksheets



In NeuronUP you can customize games and generators. Worksheets are not customizable.


1. Choose and open an activity (game or generator)



Choose a game or a generator from the NeuronUP Working Area.

Once the activity is open, select "Personalized Modes" for generators or "Personalized Stages" for games.

2. Create or edit a customization




To create a new customization, enter a name for the customization. For example, the name of the patient it will be used for.

Then click on "Continue".

To modify a customization, click on "Edit".

3. Start customizing

General features



The "general features" are the same for all games and generators.

You can choose the work mode, adapt the time, warn the patient when he/she is inactive, modify the activity instructions and add a timer.


Parameters



The "parameters" are specific to each activity.

In this case, we are in a money handling activity, so you can choose the type of currency, the amount of bills/coins and their value.


Adaptability



The "adaptability" tab allows you to choose between dragging or clicking on elements to make it easier for patients with motor impairments or those who are not good with the mouse or touch screens.

It is not available for all activities.

4. Save the customization




You can add a description to the customization.

This will make it easier for you to remember its characteristics.



To save the customization, click on "Add".

The customized activity is ready and saved




You can work directly with the customized activity.

When you want to work with it later, you will find it in the list of activity customizations.



Find your customization when you create a session

1. Go to "Sessions" > "Digital Sessions" > "New Digital Session".
2. Set up your session and add the activities by clicking on them.
3. Choose the activity for which you have created a customization. In this case *Sort the Money*.



4. Once the activity has been added, click on it on the left side of the screen, as in the following image.

5. Next, activate the "Choose where to start" filter and you will see the different customizations you have created for this activity.

6. Choose the customization you want. In this example we chose "Amelia". This way, your patient will carry out the session with the customized activity.

