

MANUAL

WORKING WITH SEVERAL PATIENTS SIMULTANEOUSLY







With NeuronUP you can get your patients to work on different devices at the same time.



To do this, **you must successively login with your password to each patient's computer** and start the session you have assigned to them.

1. Assign a session



From the menu, click on **"Sessions"** and then on **"Digital Sessions"**.

2. Choose and click on a session to assign it



Choose a session and click on it or drag it to the date you wish to assign it.

3. Choose the patients to which you are going to assign the session

neuron ¤ ,			
Piter Patient 4	ASSIGN SESSION TO PATIENT(S)		English *
🗅 Activities 🗸 🗸	9. Search Patient	Calender V	× Vew
Sessions A	NeuronUP260 Access	ch 2023 Previous Current	
Digital Sessions	10	Select day Monday The maximum duration is 15 days	Saturday
Paper Sessions	Scarlet Spinster	3/7/2023 Data - 3/7/2023	
Programs	Group 3	8 Session for home	
Extra Resources	Kristel Daft	aar aa	
	Access until Apr 1, 2023 Janelle Abbey	sa s Assigned to:	
	Alexan Carr	C Elot Callen C Janelle Abbey	
	Methe Cass		
	Eliot Callen	27 21	
	③ New Patient	Assign Session	_
← Collapse Menu			
	6	2 7	
		ELIOT CALLEN	
		Le Attention	
		Attention	

In the list of patients, choose those who will work with this session.

Note:

In order for patients to be able to work at the same time individually, it is necessary to have as many devices as patients.

For example, for 2 patients to work **at the same time**, you need 2 computers or 2 tablets.

4. Go to the first patient's device and click on "Start" in their session

From here, you will have to repeat the process on each patient's device.

			:
Not started			
Start Date * 3/1/2023	Ē	End Date *	
 Session for home 	The maximum	n duration is 15 days	
 ★ Session for home ▶ NO 	The maximun	n duration is 15 days	Save
Session for home NO Assigned by: Jacob Baker Language: English	The maximum	n duration is 15 days	Save

Once you have logged in to the first patient's computer with your user keys, click on **"Sessions"**.

Select the patient's session (in this case, Eliot Callen) and click on **"Start"**.

The patient's session starts on his or her device



From this point on, the session is secure.

The patient can only access his or her session.

Just click on **"Start"** to begin working.



Go to the second patient's device

	neuron		
	Enter		
	AUsernamo/Email* email@example.com	× •<	
	• Password ·		
	Have you forgotten your password?		
	ENTER		
	Are you still not working with NeuronUP? Free trial of NeuronUP		



Repeat the above process.

Log in with your passwords on the second patient's device.



Click on **"Digital Sessions"**.

Select the patient's session, in this case Janelle, and click on **"Start"**.

		English	•
🛱 Calendar Vie	Attention JANELLE ABBEY	×	e View
Previous Cu	Not started		
Tuesday	Start Date * End Date * 3/1/2023	Ē	Saturday
7	The maximum duration is 15 days		
, i i i i i i i i i i i i i i i i i i i	NO NO	Save	
14	 Assigned by: Jacob Baker A Language: English Time: 15 min 		
21	Start View Edit Duplicate Delete		
28	29 30	31	

Hello Janelle

Attention

Start

Start the session.

Repeat the process for all patients to whom you want to assign the session.